## Prototype Module of a Robust 18-channel Magnetometer System

J. H. Storm<sup>1</sup>, M. Burghoff<sup>1</sup>, D. Drung<sup>1</sup>, R. Körber<sup>1,2</sup>

<sup>1</sup>Physikalisch-Technische Bundesanstalt, Abbestraße 2-12, Berlin, Germany <sup>2</sup>Magnicon GmbH, Barkhausenweg 11, 22339 Hamburg , Germany

E-mail: Jan-Hendrik.Storm@ptb.de

*Abstract* — We present the prototype module of our robust multichannel SQUID magnetometer system which is intended for high-precision in biomagnetic measurements and nuclear spin precession experiments. Further demanding applications are magnetorelaxometry and ultra-low field NMR, where pulsed magnetic fields of up to 50 mT are typically applied. The system is operated in a flat bottom liquid helium dewar placed inside the Berlin Magnetically Shielded Room (BMSR-2). It contains 18 magnetometers and will be the basis for a modular multichannel system. The magnetometers are designed with Nb wire-wound flux antennas. A total of 16 small antennas (17.1 mm diameter) form a regular grid with individual channels sensitive to all three spatial directions. Two large antennas with 74.5 mm diameter sensitive in z-direction surround the grid at two different heights and allow the detection of deep sources. Each antenna is connected to the input of a thin-film Nb SQUID current sensor via a detachable contact. The SQUIDs are equipped with integrated input current limiters. Feedback into the antennas is employed to minimize crosstalk. The current sensor chip package involves a superconducting shield of Nb. The configuration of the SQUID magnetometers and the multichannel arrangement does not significantly degrade the superior magnetic environment inside the BMSR-2 which has been verified by simulations. The measured white filed noise of the small-size magnetometers was between 0.52 and 1.3 fT/VHz, and well below 1 fT/VHz for the large ones. The experimental performance of the prototype module will be discussed in detail.

*Keywords (Index Terms)* — SQUID, multichannel SQUID system, ultra-low-field NMR, MEG.